

THE SPOTV



- Designed to be placed over an existing Welch Allyn Spot Vital Sign Monitor (42NTB) - it could be an obsolete or working device.
- Allows you to push/show the following measures directly on your Welch Allyn Spot Vital signs: Systolic, diastolic, heartbeat, oxygen rate, and blood pressure.
- Compatible with low, medium, high-fidelity manikins or with a simulated/standardized patient.
- Very solid and won't damage existing devices.
- Update it via the Apple Store from the device (Spot Vital SYHM).
- Comes with a set of two rails with tape preinstalled on your SPOTV, lightning charging cable (iPhone/iPod/iPad cable) and a wall charger.
- Powered by a rechargeable Lithium-ion battery (4-hour autonomy).

THIS DEVICE IS



User-friendly



Remote-controlled via instant Bluetooth connection



Durable



Designed, manufactured and assembled in Canada



Environment-friendly
(1 device bought = 1 tree planted)



Controlled via Innov2Learn's free, easy-to-use app for all devices